



Contact Information:
Carole Lee 888-924-3800

For Immediate Release

New Mobile Dental X-ray Unveiled By Progeny, Inc.

Buffalo Grove, IL—January 2006—Progeny Dental has introduced a mobile model of their premier dental x-ray systems, the JB-70 and Preva DC.

Progeny's mobile intraoral x-ray systems provide exceptional imaging capabilities for clinical settings that require the utmost stability and maneuverability. The UL-approved systems feature heavy duty 3.5" casters with two locking wheel brakes and a design that allows the user to turn 360 degrees in place. At just 26.9" wide, Progeny Mobile Systems can easily pass through all standard doorways. A convenient handle for positioning, a corded remote control and power cord management rail make positioning and storage simple and efficient.

Preva DC Mobile combines the flexibility of adjustable kV, mA and time techniques with a 0.4mm focal spot, making it the ideal choice for digitally-based practices.

JB-70 Mobile features a 0.7mm focal spot, 70kV and 7mA technique factors and advanced AC technology optimized to provide the highest quality images using both film and digital sensors.

For additional information, visit www.progenydental.com or call 888-924-3800, Ext. 254.

-END-