

Operating System Requirements For Progeny VisionDX Sensor System

- **Vista – 32 bit**
- **Windows XP™ Pro with Service Pack 2 – 32 bit**
- **Apple MacBook Pro with Microsoft Windows Platform
(Contact Progeny, Inc. for complete details)**
- **Hard Drive: 60 GB minimum (2GB RAM free for installation)**
- **Processor speed 2GHz or better**
- **RAM: 2GB free disk space minimum**
- **CD ROM Drive: 48X speed**
- **Video display adapter: 512K, 256K acceptable**
- **Display: 1024 x 768**
- **32 bit, true color recommended**
- **USB 2.0 port available**
- **Back up device recommended (DVD, External Hard Drive,
other)**

Updated 10/23/2008

Subject to change without notice.

Contact Progeny, Inc. 888-924-3800, Ext. 2