



Progeny
A Midmark Company
675 Heathrow Drive
Lincolnshire, Illinois 60069-4206
1-847-415-9800
Fax 847-415-9801
progenydenal.com

Progeny Imaging Software System Requirements

Operating System	Microsoft Windows 7 [®] Ultimate or Professional (32 & 64 bit) or Microsoft Vista [®] Business or Ultimate (32 & 64 bit) or Microsoft Windows XP [®] Professional (with SP2 or SP3) or Intel [®] -based Apple [®] running Windows Platform.
Processor	Intel [®] (or x86 compatible) Single Core 2 GHz or better (or Dual Core at 1.7GHz recommended)
Memory	2 GB RAM (minimum), 3 GB RAM (recommended)
Storage	80GB Hard Drive (minimum)
CD/DVD	48x CDRW (minimum), 16x DVD+/-RW (recommended)
Video	32 bit, 1024 X 768 Resolution Capable (minimum)
Display	1024 X 768, 32 true bit color (recommended)
Ports	Two USB 2.0

A back-up storage device is recommended (DVD, External Hard Drive, other)

Note: Microsoft Windows 7[®] Home Premium, Microsoft Windows Vista[®] Home Basic or Home Premium and Microsoft XP[®] Home and Tablet PC Editions are NOT recommended.

For additional information, contact Progeny Technical Support at 847-415-9800 option 2.

U.S. Patent No. 5.179.579

Effective December 8, 2009. Subject to change without notice.

ML-00009 Rev. B